



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Carbolite-Gero® TSO3 11/600 1100degrees, 600mm EPC3016P1 PID controller TSO3 11/600 1100degrees, 600mm

**€18,010.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

TSO3 11/600 1100°C, 600mm Heated length, three zone rotary reactor tube furnace (220 - 240 V three phase) (with setpoint adjustment) CC-T1 Program Controller and Data Recorder with USB and Ethernet - Single zone control - 4.3" touch screen - 10 programs - 24 segments per program (e.g. 12 ramps and 12 dwells) - Max. 2 control tracks (e.g. for optional solenoid valves) - USB port on the front of the control box - Ethernet port on the back of the control box - Data recording in csv format (comma separated variables) on customer USB stick - Real time clock with calendar - Date and time for program start programmable - Program status display - Output level display of heaters - Language selection: English, German, French, Italian, Spanish, Chinese - Password protected user levels: 1. Operator can set temperature setpoints / run preset programs / start data logging 2. Supervisor can additionally change and create programs / acknowledge alarms 3. Administrator can additionally change language settings, time and date / change data logging presets / change calibration settings and communication settings Standard TSO/600 quartz tank with smooth inner surface Optional overtemperature protection digital (00053325) Upgrade TSO/600 to quartz tank with fluted inner surface

## Additional Information

---

No.	CX-0407
Manufacturer (Brand)	Carbolite-Gero
Transport temperature	Room temperature

