



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Fritsch® rotor and holder for sieve trays with steel knives

€2,871.00
plus VAT &
Shipping

Product Images

Description

Cutting rotor made of stainless steel - for Variable Speed Rotor Mill PULVERISETTE 14 premium line consisting of cutting rotor with cooling ribs and sieve pan support made of stainless steel and rotor blades and counter blades made of hardened stainless steel. The stainless steel cutting rotor is ideal for all standard tasks of pre-shredding and fine shredding by cutting and shearing as well as for fast problem-free shredding of fibrous materials and plastics. Choose 316L stainless steel screen bowls with trapezoidal or round holes for the desired final fineness.

Additional Information

No.	FR-0954
Manufacturer (Brand)	Fritsch
Transport temperature	Room temperature

