



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Fritsch® set = 4 pieces sieve trays 1,5 mm round hole P-14 pl stainless steel - cutting mill**

**€134.00  
plus VAT &  
Shipping**

**Product Images**

---

## Description

---

Set of sieve bowls 1.5 mm round hole made of stainless steel In just a few steps, the Variable Speed Rotor Mill PULVERISETTE 14 premium line becomes a cutting mill for the quick, problem-free pre-shredding of fibrous materials and plastics. Together with the cutting rotor with cooling ribs, simply insert the suitable accessories such as, labyrinth disk, collecting vessel, holder for sieve trays and a pair of sieve trays - the comminution then takes place by cutting and shearing. The desired final fineness is also determined by the selected sieve tray with trapezoidal or round perforation. Sieve trays with round perforations are recommended for shredding brittle material and medium fineness with a narrow grain width band. Sieve trays with trapezoidal perforations are also available. These are suitable for rapid size reduction of medium-hard to soft materials; the trapezoidal perforation produces an additional shearing effect. Sieve trays are also available with other perforations.

## Additional Information

---

No.	FR-0118
Manufacturer (Brand)	Fritsch
Transport temperature	Room temperature

