



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Fritsch® dividing head cpl. 3 samples ALU head PTFE coated

€4,156.00
plus VAT &
Shipping

Product Images

Description

PTFE-coated aluminum is recommended for the division of aggressive solids and suspensions that attack plastic POM and aluminum. Dividing ratio: 1:30 Number of dividing quantities: 3 max. feed size: 2.5 mm max. sample quantity: 300 ml Division: dry/wet Useful volume of sample vials: 15, 20, 30 ml Equipment: This dividing head is supplied as standard with 3 sample vials 15 ml with screw cap and funnel 10 mm Ø as well as 1 funnel 5 mm Ø each for dry and wet division (included in the price of the dividing head). This dividing head is also available in plastic POM, article number 27.6150.00 for dividing materials and suspensions that do not produce abrasion on the surface.

Additional Information

No.	FR-0207
Manufacturer (Brand)	Fritsch
VGKL number	549827215
Transport temperature	Room temperature

