



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

**Fritsch® conversion kit for shredding larger quantities Collecting vessel with gr. recept. - from 01**

**€1,427.00  
plus VAT &  
Shipping**

**Product Images**

---

## Description

---

consisting of pan with outlet and flanged nylon support bag (outside) for paper filter bag (inside). Regardless of the specific configuration with rotor, sieve ring, beater bar or pin insert, the PULVERISETTE 14 can grind larger quantities with this special conversion kit: Conversion kit consisting of a special collecting vessel with outlet and a flanged-on nylon support bag for easily replaceable paper filter bags. This allows sample quantities of up to 1 l to be ground in one operation without having to open the grinding chamber and empty the collecting pan in between. Our tip : When grinding very light materials, the material discharge can be significantly improved by the high air circulation with the conversion kit for grinding larger quantities - even if you are grinding smaller quantities! Another tip : When grinding temperature-sensitive samples, the conversion kit for grinding larger quantities with its large nylon support bag ensures a higher air throughput and thus better cooling.

## Additional Information

---

No.	FR-0067
Manufacturer (Brand)	Fritsch
VGKL number	109801000
Transport temperature	Room temperature

