



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



KERN® load cell, 50kg (~500N) "s",

Price on
request

plus VAT &
Shipping

Product Images



Description

--Accuracy class according to OIML C3 --RoHS compliant --Dust and splash protection IP68 --Rust-free steel --
Application area: mass and force measurement --Suitable for hanging scales, silo scales, force test benches and
other various scales --Nominal rating: 2.0 ± 0.004 --CS P2 50–250 kg miniature S load cells/S load cells made of
stainless steel

Additional Information

No.	KP-9644
Manufacturer (Brand)	KERN
Transport temperature	Room temperature

