



neoLab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## KERN® load cell, bending bar load cell 250kg

**Price on  
request  
plus VAT &  
Shipping**

### Product Images



## Description

---

Accuracy according to OIML R60 C3; CE and RoHS compliant; Dust and splash water protection IP67; (according to EN 60529); Nickel-plated steel; Field of application: Mass as well as compressive force measurements under harsh environmental conditions; Suitable for motor vehicle scales, hopper scales, motor vehicle test equipment, test benches; 4-wire connection; Nominal sensitivity: 3 mV/V;;

## Additional Information

---

No.	KP-8839
Manufacturer (Brand)	KERN
Transport temperature	Room temperature

