



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lauda® Proline PV 36 see-through therm.
208-220 V, 60 Hz**

**€12,400.00
plus VAT &
Shipping**

Product Images

Description

LAUDA Proline PV 36 see-through therm. 208-220 V; 60 Hz Performance features: Bath/circulator thermostat with state-of-the-art microprocessor technology and built-in cooling coil Easy-to-read green LED display for temperature User-friendly menu navigation and simplest 3-button operation EasyUse system for easiest handling and operation of the entire unit SelfCheck assistant for system diagnostics Fully electronic continuous controller with PID behavior for internal and external control PowerAdapt system for optimally adapted max. Heating power without power supply overload Low level protection and adjustable overtemperature protection with acoustic alarm for operation with flammable and non-flammable liquids LAUDA Varioflex pump (pressure pump) with 8 selectable power levels Retrofittable with up to 2 additional interface modules (RS 232/485, Profibus, analog or contact modules, Ethernet-USB module) Measuring chamber with large viewing windows at the front and rear Pump connections at the rear and side, integrated bypass Bath vessel made of stainless steel (insulated, with carrying handles and drain cock) Technical data (according to DIN 12876): Working temperature range: 30 ... 230 °C Working temperature range with water cooling: 20 ... 230 °C Operating temperature range: 0 ... 230 °C Ambient temperature range: 5 ... 40 °C Temperature stability: 0.01 ± K Heating power max.: 3.3 kW Current consumption: 16 A Power consumption max.: 3.5 kW Pump flow max. (pressure): 25 L/min Bath volume min. / max.: 28.0 / 36.0 L Bath size (WxDxH): 585 x 135 x 320 mm Dimensions (WxDxH): 1040 x 282 x 590 mm Weight: 47 kg Power supply: 208-220 V; 60 Hz Power plug: Power cord with plug (NEMA 6-20P) Standard accessories: 4 union nuts, 4 sealing plugs 2 hose olives 13 mm for pump connection 2 hose olives for cooling coil

Additional Information

No.	LA-2439
Manufacturer (Brand)	Lauda
Transport temperature	Room temperature

