



# Liebherr® SmartCoolingHub

## Product Images

## Description

The SmartCoolingHub is the basis for networking your professional Liebherr appliances. Up to 20 networkable Liebherr appliances can be connected by cable. The SmartCoolingHub is connected to your network and sends the device data to the Liebherr cloud, which is clearly displayed in the SmartMonitoring dashboard. The SmartCoolingHub also has a buffer memory and records all device data for up to 48 hours if the Internet connection is interrupted. Connection WLAN 2.4 GHz or LAN 100Mbit Up to 20 Liebherr devices can be networked with serial interface RS 485 Up to 48 hours of data storage in case of interruption of the Internet connection Space-saving design, wall mounting possible Status LED shows the operating status

## Additional Information

No.	LH-0195
Manufacturer (Brand)	Liebherr

