



## Lovibond® Fluoride Standard

### Product Images

---



## Description

Calibration Standard Solution 1 mg/l, 30 ml For use with SPADNS Reagent Reagent solutions usually do not consist of one preparation, but of several components that must be added to the sample one after the other. Since both size and number of drops have a decisive influence on the color complex to be formed, the dosage must be extremely accurate. The shelf life of liquid reagents is limited by temporary contact with atmospheric oxygen when the bottles are opened and by unsuitable storage conditions (exposure to sunlight, higher temperatures).

## Additional Information

No.	CV-0003
Manufacturer (Brand)	Lovibond

