



Lovibond® Hydrazine Test Powder

Product Images



Description

Powder, 30 g The simple and fast form of administration make the powder reagents popular reagents for water analysis in many countries. The Lovibond ® powder reagents program offers familiar users an alternative for existing measuring systems. The powder reagents are manufactured under the same qualitative requirements that have made Tintometer GmbH, Lovibond ® Water Testing so successful in the reagent tablet sector for decades. Parameters from aluminum to chlorine to sulfate, are known worldwide parameters that can be covered with the powder reagents program.

Additional Information

No.	CV-0191
Manufacturer (Brand)	Lovibond

