



Lovibond® Reagent Tablet Copper/Zinc LR

Product Images



Description

Reagent tablet Copper/Zinc LR for the determination of zinc. VE=100 tablets. (VE 250 tablets on request) With the right tablet for fast results Water analysis could not be simpler. Precise dosing, enormous shelf life and easy handling make the reagent tablet the most popular indicator system and the most economical way to measure a parameter. The reagent tablet defies almost all climatic conditions and external influences in its aluminum blister packaging. In addition, there is no safety risk for the user in the application.

Additional Information

No.	CV-0414
Manufacturer (Brand)	Lovibond

