



## Lovibond® reagent tablet CyA Test

### Product Images

---



## Description

Reagent tablet CyA Test Test-Parameter Stabilzier (cyanuric acid) in blister pack. VE=100 tablets. (VE 250 tablets on request) With the right tablet for fast results Water analysis could not be easier. Precise dosing, enormous shelf life and easy handling make the reagent tablet the most popular indicator system and the most economical way to measure a parameter. The reagent tablet defies almost all climatic conditions and external influences in its aluminum blister packaging. In addition, there is no safety risk for the user in the application.

## Additional Information

|                      |          |
|----------------------|----------|
| No.                  | CV-0405  |
| Manufacturer (Brand) | Lovibond |

