



Lovibond® reagent tablet DPD No.4

Product Images

Description

Reagent tablet DPD No.4, test parameter active oxygen in blister pack. VE=100 tablets (VE 250 and 500 tablets on request) With the right tablet for fast results Water analysis could not be easier. Precise dosing, enormous shelf life and easy handling make the reagent tablet the most popular indicator system and the most economical way to measure a parameter. The reagent tablet defies almost all climatic conditions and external influences in its aluminum blister packaging. In addition, there is no safety risk for the user in the application.

Additional Information

No.	CV-0404
Manufacturer (Brand)	Lovibond

