



Lovibond® reagent tablet EDTA

Product Images



Description

Reagent tablet EDTA for zinc. VE=100 tablets (VE 250 tablets on request) With the suitable tablet to the fast result Water analysis cannot be simpler. Precise dosing, enormous shelf life and easy handling make the reagent tablet the most popular indicator system and the most economical way to measure a parameter. The reagent tablet defies almost all climatic conditions and external influences in its aluminum blister packaging. In addition, there is no safety risk for the user in the application.

Additional Information

No.	CV-0411
Manufacturer (Brand)	Lovibond

