



## Lovibond<sup>®</sup> reagent tablet Silica PR

### Product Images

---



## Description

---

Reagent tablet Silica PR for silicon dioxide (in the presence of phosphate), PU = 100 tablets (PU 250 tablets on request) Successor of: 4513150BT With the right tablet for a fast result Water analysis cannot be easier. Precise dosing, enormous shelf life and easy handling make the reagent tablet the most popular indicator system and the most economical way to measure a parameter. The reagent tablet defies almost all climatic conditions and external influences in its aluminum blister packaging. In addition, there is no safety risk for the user in the application.

## Additional Information

---

No.	CV-0422
Manufacturer (Brand)	Lovibond

