



Lovibond® Reagent Tablet Sulfites LR

Product Images

Description

Reagent tablet Sulfite LR for sulfite PU = 100 tablets With the right tablet for a fast result Water analysis cannot be simpler. Precise dosing, enormous shelf life and easy handling make the reagent tablet the most popular indicator system and the most economical way to measure a parameter. The reagent tablet defies almost all climatic conditions and external influences in its aluminum blister packaging. In addition, there is no safety risk for the user in the application.

Additional Information

No.	CV-0447
Manufacturer (Brand)	Lovibond

