



Lovibond® VARIO Sulfa 4 F10

Product Images

Description

VARIO Sulphate 4 F10, Powder Packs, PU = 100 Powder Packs The simple and fast form of administration make the powder reagents popular reagents for water analysis in many countries. The Lovibond ® Powder Reagents range offers familiar users an alternative for existing measuring systems. The powder reagents are manufactured under the same qualitative requirements that have made Tintometer GmbH, Lovibond ® Water Testing so successful in the reagent tablet sector for decades. Parameters from aluminum to chlorine to sulfate, are known worldwide parameters that can be covered with the powder reagents program.

Additional Information

No.	CV-0457
Manufacturer (Brand)	Lovibond

