



Lovibond® Vario Phosphate Reagent F10

Product Images



Description

VARIO Phosphate Reagent F10, Powder Pack, PU = 100 Powder Packs The simple and fast form of administration make the powder reagents popular reagents for water analysis in many countries. The Lovibond Powder Reagents range offers familiar users an alternative for existing measuring systems. The powder reagents are manufactured under the same qualitative requirements that have made Tintometer GmbH, Lovibond Water Testing so successful in the reagent tablet sector for decades. Parameters from aluminum to chlorine and sulfate are known worldwide and can be covered with the powder reagents program.

Additional Information

No.	CV-0456
Manufacturer (Brand)	Lovibond

