




Macherey & Nagel® Chromab.S.SiOH-H₂SO₄/SA,3mL,500/500mgBIG

Product Images



Description

CHROMABOND columns SiOH-H2SO4/SA BIGpack Volume: 3 mL, filling quantity: 500/500 mg Material: PP, with PE filter elements packed in 5 bags à 50 pcs. Packing à 250 pcs. EQ Exempted quantity or __UN 3316 Chemistry test kit 9 II 0,125 KG ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 KG CAO 960 10 KG

Additional Information

No.	ML-8223
Manufacturer (Brand)	Macherey & Nagel

