



## Macherey & Nagel® NucleoTraPCR (100)

### Product Images

---



## Description

NucleoTraP CR (100) 100 preps for PCR clean-up - NucleoTraP CR Suspension, buffers pack of 100 \_\_UN 3211 Perchlorates, inorganic, aqueous solution, n.o.s. 5.1 III 0.100 L ADR/GGVSEB 01 (E) LQ PAX 551 2.5 L CAO 555 30 L Q: 0.040 / 0.004

## Additional Information

|                      |                  |
|----------------------|------------------|
| No.                  | ML-4195          |
| Manufacturer (Brand) | Macherey & Nagel |

