



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Macherey & Nagel® VISOCOLOR HE Cyanide High Sensitivity Test Kit - Refill - Measuring range: 0.002-0.040 mg/L CN-

**€69.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

VISOCOLOR HE Cyanide High Sensitivity Test Kit - Refill - Measuring range: 0.002-0.040 mg/L CN- Pg. à 50 determinations \_\_UN 3316 Chemistry Test Kit 9 II 0.040 L ADR/GGVSEB M11 (E) ADR 3.3.1/251: LQ PAX 960 10 Kg CAO 960 10 Kg

## Additional Information

---

No.	ML-6680
Manufacturer (Brand)	Macherey & Nagel
Transport temperature	Room temperature

