



Rötzmeier® Watering Tank 4l: TB4

Product Images



Description

Impregnation container 4l: TB4 Safety container made of stainless steel, made in Germany (Salzkotten) Container made of high quality and durable stainless steel - wall thickness 0.8 mm For controlled impregnation of rags, e.g. impregnation with solvents for cleaning parts Minimizes consumption High stability of the container Flame arrestor in the impregnation plate Impregnation plate spring-loaded Laboratory pot, stainless steel, vessel, vessels, safety stand vessels, container, stainless steel container

Additional Information

No.	RZ-0185
Manufacturer (Brand)	Rötzmeier

