



Sartorius® Microsart Funnel 250

Product Images



Description

Microsart® Funnels are sterile plastic funnels. They allow rapid performance of filtration in routine testing of water, beverages and food, pharmaceutical and cosmetic products. After placing a Sartorius Stedim Biotech grid membrane on a sterile stainless steel filter support, the Microsart® Funnel is easily and conveniently placed and the sample is vacuum filtered. The funnel is made of polypropylene and is therefore elastic enough to ensure an optimal seal by means of a click-fit closure. A graduation with markings allows sufficiently accurate volume dosing. The large inner diameter allows a high flow rate and the optimized inlet of the inner wall allows a perfect rinsing of the system. No residual fluid remains in the filter funnel. Reliable results - a new sterile filter funnel for each test safely avoids secondary contamination Time saving - no more wasted time due to flaming Easy handling - no more handling of hot filter funnels

Additional Information

No.	ST-3562
Manufacturer (Brand)	Sartorius

