



## Sartorius® Vivapure AdenoPACK™ 500

### Product Images

---



### Description

---

Vivapure AdenoPACK 500 is the premium version of Adeno PACK 100 for virus purification. In only 2 hours,  $3 \times 10^{13}$  adenovirus particles are purified and concentrated from 500 ml cell culture. This ready-to-use kit is easy and convenient to use with a laboratory pump, ensuring optimal flow control with minimal time consumption. This easy-to-use product replaces long and inefficient cesium chloride density measurement methods. Pump-operated kit Easy and controlled purification of adenoviruses from large cell culture volumes Virus purification based on membrane adsorber technology High-speed process ready in 2 hours Convenient, easy-to-use alternative to CsCl Density Gradients Suitable to handle 10 to 25 plates (up to 500 ml) Flexible, most cost-effective kit on the market for purification of  $1-3 \times 10^{13}$  virus particles

# Additional Information

|                      |           |
|----------------------|-----------|
| No.                  | ST-6566   |
| Manufacturer (Brand) | Sartorius |

