



Nalgene® stopcocks

Product Images



Description

Nalgene stopcocks Localized liquid transfer systems for accurate filling or dispensing. These stopcocks provide reliable and accurate flow control of liquids. Type: Stopcock Plug Opening: 2 mm Quantity per package: 1 Quantity per carton: 6 For use with (devices): In-House Fluid Transfer System Plug Material: Polypropylene Unit Size: Case of 6

Additional Information

No.	CD-1776
Manufacturer (Brand)	Nalgene
VGKL number	278764601

