



## Nalgene® Acrylic Parafilm Dispenser

### Product Images

---



16

## Description

Nalgene Acrylic Parafilm Dispenser The Thermo Scientific Nalgene Acrylic Parafilm Dispenser Box is designed to securely hold Parafilm rolls for easy removal Quantity per box: 1 Quantity per package: 1 Height (metric): 152 mm Length (metric): 152 mm Width (metric): 121 mm Description: acrylic parafilm dispenser box Item description: AD, 152 times; 121 times; 152 mm Unit Size: Each

## Additional Information

|                      |         |
|----------------------|---------|
| No.                  | CD-1612 |
| Manufacturer (Brand) | Nalgene |

