



Nalgene® Capitol Vial Collection Cup

Product Images

Description

Capitol Vial Collection Beakers Thermo Scientific Capitol Vial Collection Beakers are suitable for conventional and on-site testing.
Material: Plastic Type: Capitol Vial Collection Cup Volume (metric): 236 mL Description: collection beaker with temperature strip; 8 oz. (236 mL) Contains: Temperature Strips Quantity per Case: 500 Volume (Anglo): 8 oz. Unit Size: Case of 500

Additional Information

No.	CD-2737
Manufacturer (Brand)	Nalgene

