



Nalgene® Capitol Vial Collection Cup

Product Images

Description

Capitol Vial Collection Beakers Thermo Scientific Capitol Vial Collection Beakers are suitable for conventional and on-site testing.
Material: Plastic Type: Capitol Vial Collection Cup Volume (metric): 103.5 mL Description: collection beakers with temperature strips;
3.5 oz. (103.5 mL) Contains: Temperature strips Quantity per box: 350 Volume (anglo): 3.5 oz. Unit Size: Case of 350

Additional Information

No.	CD-2738
Manufacturer (Brand)	Nalgene

