



Nalgene® Capitol Vial Collection Cup

Product Images

Description

Capitol Vial Collection Cup Thermo Scientific Capitol Vial Collection Cups are suitable for conventional and on-site testing. Type: Specimen Collection Cup Volume (metric): 236 mL Description: Beaker, 8 oz., with Celsius Temperature Strip Contains: Temperature Strip in degrees Celcius Number per Case: 500 Volume (anglo): 8 oz. Unit Size: Case of 500

Additional Information

No.	CD-2733
Manufacturer (Brand)	Nalgene

