



Nalgene® Closed-Dome Polypropylene Tanks

Product Images



Description

Nalgene Closed-Dome Polypropylene Tanks Storage, mixing and dispensing of reagents in Thermo Scientific™ Nalgene™ closed dome polypropylene tanks. Certification/Compliance: Comply with FDA Reg. 177.1520 and USP Class VI Outer Diameter (Anglo-American): 18.5 in. Type: Closed dome tank Capacity (anglo): 30 gal. Capacity (metric): 115 L Height (Anglo): 38.625 in. Unit Size: Each

Additional Information

No.	CD-0678
Manufacturer (Brand)	Nalgene

