



Nalgene® Closed-Dome Polypropylene Tanks

Product Images



Description

Nalgene Closed-Dome Polypropylene Tanks Storage, mixing and dispensing of reagents in Thermo Scientific™ Nalgene™ closed dome polypropylene tanks. Certification/Compliance: Comply with FDA Regulation 177.1520 and USP Class VI Outer Diameter (Anglo-American): 18.5 in. Outer Diameter (metric): 47 in. Type: Closed dome container Capacity (Anglo-American): 20 gal. Capacity (metric): 75 l Height (Anglo-American): 35 in. Height (metric): 81 cm Unit Size: Each

Additional Information

No.	CD-0677
Manufacturer (Brand)	Nalgene

