



Nalgene® Closed-Dome Polypropylene Tanks

Product Images



16

Description

Nalgene Closed-Dome Polypropylene Tanks Storage, mixing and dispensing of reagents in Thermo Scientific™ Nalgene™ closed dome polypropylene tanks. Certification/Compliance: Comply with FDA Reg. 177.1520 and USP Class VI Outer Diameter (Anglo-American): 22 in. Type: Closed dome tank Capacity (Anglo): 55 gal. Capacity (metric): 210 L Height (Anglo): 39.5 in. Unit Size: Each

Additional Information

No.	CD-0679
Manufacturer (Brand)	Nalgene

