



Nalgene® Depyrogenated Sterile Empty Vials

Product Images



Description

Depyrogenated Sterile Empty Vials Thermo Scientific depyrogenated sterile empty vials are available in sizes from 1 to 100 ml.
Description: Clear glass tube, assembled with stopper and aluminum flare seal Outer diameter (anglo): 0.5 in. Outer diameter (metric): 12.7 mm Material: type 1, borosilicate glass Volume (metric): 2 ml Height (anglo): 1.313 in. Quantity per case: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Closure Material: aluminum Closure Type: flare cap (assembled) with stopper Quantity: 100 boxes Certification/Compliance: certificate of sterility and pyrogen test are included Type: vials Closure Size: 13 mm For use with (devices): Aseptic protocol Unit Size: Each

Additional Information

No.	CD-3707
Manufacturer (Brand)	Nalgene

