



Nalgene® Depyrogenated Sterile Empty Vials

Product Images



Description

Depyrogenated Sterile Empty Vials Thermo Scientific depyrogenated sterile empty vials are available in sizes from 1 to 100 ml.
Description: Amber glass vial pre-assembled with stopper and aluminum crimp seal Outer diameter (anglo): 0.813 in. Volume (metric): 10 mL Height (anglo): 2 in. Quantity per box: 50 Volume (anglo): 0.34 oz. Total Volume (metric): 10 mL Closure Material: Aluminum Closure Type: Crimp Cap (Assembled) with Stopper Quantity: 50 Cases Certification/Compliance: Certificates of Sterility and Pyrogen Test included Type: Vial Closure Size: 20 mm For Use with (Devices): Aseptic Protocol Unit Size: Each

Additional Information

No.	CD-3705
Manufacturer (Brand)	Nalgene

