



Nalgene® Depyrogenated Sterile Empty Vials

Product Images



Description

Depyrogenated Sterile Empty Vials Thermo Scientific depyrogenated sterile empty vials are available in sizes from 1 to 100 ml.
Description: Amber glass vial pre-assembled with stopper and aluminum crimp seal Outer diameter (Anglo-American): 1.313 in.
Volume (metric): 30 mL Height (anglo): 2.5 in. Quantity per box: 50 Volume (anglo): 1.014 oz. Total Volume (metric): 30 mL Closure
Material: Aluminum Closure Type: Crimp Cap (Assembled) with Stopper Quantity: 50 Cases Certification/Compliance: Certificates of
Sterility and Pyrogen Test included Type: Vial Closure Size: 20 mm For Use with (Devices): Aseptic Protocol Unit Size: Each

Additional Information

No.	CD-3710
Manufacturer (Brand)	Nalgene

