



Nalgene® Depyrogenated Sterile Empty Vials

Product Images



Description

Depyrogenated Sterile Empty Vials Thermo Scientific depyrogenated sterile empty vials are available in sizes from 1 to 100 ml.
Description: Clear glass vials, assembled with stopper and aluminum flare seal Outer diameter (Anglo-American): 1.5 in. Volume (metric): 50 ml Height (Anglo): 2.688 in. Quantity per carton: 50 Volume (anglo): 1.7 oz. Total Volume (metric): 50 mL Closure Material: aluminum Closure Type: flare cap (assembled) with stopper Quantity: 50 cartons Certification/Compliance: certificate of sterility and pyrogen test enclosed Type: vials Closure Size: 20 mm For use with (devices): Aseptic protocol Unit Size: Each

Additional Information

No.	CD-3711
Manufacturer (Brand)	Nalgene

