



Nalgene® Depyrogenated Sterile Empty Vials

Product Images



Description

Depyrogenated Sterile Empty Vials Thermo Scientific depyrogenated sterile empty vials are available in sizes from 1 to 100 ml.
Description: Amber glass vial pre-assembled with stopper and aluminum crimp seal Outer diameter (anglo): 0.5 in. Volume (metric): 2 mL Height (anglo): 1.313 in. Quantity per carton: 100 Volume (anglo): 0.068 oz. Total Volume (metric): 2 mL Closure Material: Aluminum Closure Type: Crimp Cap (Assembled) with Stopper Quantity: 100 Cases Certification/Compliance: Certificates of Sterility and Pyrogen Test included Type: Vial Closure Size: 13 mm For Use with (Devices): Aseptic Protocol Unit Size: Each

Additional Information

No.	CD-3708
Manufacturer (Brand)	Nalgene

