



Nalgene® Depyrogenated Sterile Empty Vials

Product Images



Description

Depyrogenated Sterile Empty Vials Thermo Scientific depyrogenated sterile empty vials are available in sizes from 1 to 100 ml.
Description: Clear glass vial pre-assembled with stopper and aluminum crimp seal Outer diameter (angle): 1.44 in. Volume (metric): 30 ml Height (angle): 2.72 in. Quantity per carton: 50 Volume (angle): 1.014 oz. Total Volume (metric): 30 mL Closure Material: aluminum Closure Type: flare cap (assembled) with stopper Quantity: 50 cartons Certification/Compliance: certificate of sterility and pyrogen test enclosed Type: vials Closure Size: 20 mm For use with (devices): Aseptic protocol Unit Size: Each

Additional Information

No.	CD-3709
Manufacturer (Brand)	Nalgene

