



Nalgene® Depyrogenated Sterile Empty Vials

Product Images



Description

Depyrogenated Sterile Empty Vials Thermo Scientific depyrogenated sterile empty vials are available in sizes from 1 to 100 mL. Description: Clear glass vial pre-assembled with stopper and aluminum crimp seal Outer diameter (Anglo-American): 1.75 in. Volume (metric): 100 mL Height (anglo): 3.75 in. Quantity per carton: 50 Volume (anglo): 3.4 oz. Total Volume (metric): 100 mL Closure Material: aluminum Closure Type: flare cap (assembled) with stopper Quantity: 50 cartons Certification/Compliance: certificate of sterility and pyrogen test enclosed Type: vials Closure Size: 20 mm For use with (devices): Aseptic Protocol Unit Size: Each

Additional Information

No.	CD-3702
Manufacturer (Brand)	Nalgene

