



Nalgene® Depyrogenated Sterile Empty Vials

Product Images



Description

Depyrogenated Sterile Empty Vials Thermo Scientific depyrogenated sterile empty vials are available in sizes from 1 to 100 ml.
Description: Clear glass vial pre-assembled with stopper and aluminum crimp seal Outer diameter (anglo): 0.781 in. Volume (metric): 5 ml Height (anglo): 1.5 in. Quantity per box: 50 Volume (anglo): 0.17 oz. Total Volume (metric): 5 mL Closure Material: aluminum Closure Type: flare cap (assembled) with stopper Quantity: 50 cartons Certification/Compliance: certificate of sterility and pyrogen test enclosed Type: vials Closure Size: 20 mm For use with (devices): Aseptic protocol Unit Size: Each

Additional Information

No.	CD-3714
Manufacturer (Brand)	Nalgene

