



Nalgene® EFT quick-acting taps

Product Images



Description

Nalgene EFT Quick-Action Spigots Can be used with: DS2327 Series Fluorinated Balloon Bottles Description: Quick-Action Spigot ETFE Fits 1 -1/8 Inch Tap Attachment for Fluorinated Balloon Bottles 1 -1 /8 Inch Inner Diameter (Anglo) Tubing: 5/8 Inch Inner Diameter (Metric) Tubing: 16 mm Type: Spigot Quantity Per Case: 2 Quantity Per Package: 1 For Use With (Devices): DS2327 series fluorinated balloon cylinders Unit Size: Case of 2

Additional Information

No.	CD-1775
Manufacturer (Brand)	Nalgene
VGKL number	228721622

