



Nalgene® HDPE Cylinder Round Bottles

Product Images

Description

HDPE Cylinder Round Bottles Collect water samples and other low-viscosity liquids to be analyzed for inorganic trace analytes with the original EPA-specified bottle for liquid sampling. Capacity: 16 oz., 500 mL Capacity (metric): 500 mL Capacity (anglo): 16 oz. Certification Parameters: None Description: HDPE Cylinder Round Quantity per case: 120 Packaging: Case of 120 Replacement closure: 067-2810WBP Transparency: Translucent Certified: None Closure Color: White Outer Diameter (Anglo): 2.44 in. Fill Level: Unprocessed Volume to Rim (Anglo): 17.3 oz. Volume to Rim (metric): 511 mL Certification/Compliance: None Closure Liner: Foam/LDPE Closure Size: 28-410 mm Closure Material: Polypropylene Height (Anglo): 7.87 in. Shape: Round Unit Size: Case of 120

Additional Information

No.	CD-0056
Manufacturer (Brand)	Nalgene

