



Nalgene® HDPE Cylinder Round Bottles

Product Images

Description

HDPE Cylinder Round Bottles Collect water samples and other low-viscosity liquids to be analyzed for inorganic trace analytes with the original EPA-specified bottle for liquid sampling. Capacity: 4 oz., 125 mL capacity (metric): 125 mL capacity (anglo): 4 oz. Certification parameters: Metals and water quality, Cyanide, Fluoride Description: HDPE Cylinder Round Number per box: 48 Packaging: Case of 48 Replacement cap: 067-2410WBP Transparency: Translucent Certified: Metals and water quality, CN, F Closure Color: White Outer Diameter (Anglo): 1.65 in. Fill Level: Certified Volume to Rim (Anglo): 4.4 oz. Volume to Rim (metric): 131 mL Certification/Compliance: Certified containers meet or exceed EPA performance-based standards for metals, cyanide, and fluoride; manufactured under an ISO certified quality management system. Closure Liner: Foam/LDPE Closure Size: 24-410 mm Closure Material: Polypropylene Height (Anglo): 4.68 in. Shape: Round Unit Size: Case of 48

Additional Information

No.	CD-0955
Manufacturer (Brand)	Nalgene

