



Nalgene® HDPE Cylinder Round Bottles

Product Images

Description

HDPE Cylinder Round Bottles Collect water samples and other low-viscosity liquids to be analyzed for inorganic trace analytes with the original EPA-specified bottle for liquid sampling. Capacity: 32 oz., 1000 mL capacity (metric): 1000 mL capacity (anglo): 32 oz. Certification parameters: Metals and water quality, Cyanide, Fluoride Description: HDPE Cylinder Round Inner diameter (anglo): 0.81 in. Quantity per box: 12 Packaging: Case of 12 Replacement cap: 067-2810WBP Transparency: Translucent Certified: Metals and water quality, CN, F Closure Color: White Outer Diameter (Anglo): 3.11 in. Fill Level: Certified Volume to Rim (Anglo): 34.8 oz. Volume to Rim (metric): 1028 mL Certification/Compliance: Certified containers meet or exceed EPA performance-based standards for metals, cyanide, and fluoride; manufactured under an ISO certified quality management system. Closure Liner: Foam/LDPE Closure Size: 28-410 mm Closure Material: Polypropylene Height (Anglo): 9.48 in. Shape: Round Unit Size: Case of 12

Additional Information

No.	CD-0958
Manufacturer (Brand)	Nalgene

