



Nalgene® Nalgene

Product Images



Description

Nalgene Thermo Scientific Array: 4 x 4 Capacity: 80 low-profile , 48 standard , 16 deep well Shape: Rectangle Type: Microplate Rack
Quantity per box: 1 Quantity per package: 1 Dimensions (L x W x H): 15 x 5.5 x 8.5 in. (38.7 x 13.9 x 21.6cm) Length (Anglo-American): 15 in. Width (anglo): 5.5 in. Height (anglo): 8.5 in. Unit Size: Each

Additional Information

No.	CD-1541
Manufacturer (Brand)	Nalgene

