



Nalgene® PETG Biotainer Bottles

Product Images



15

Description

Nalgene PETG Biotainer Bottles Save on cleaning, sterilization and related validation costs compared to reusable bottles. Capacity (metric): 1 L Capacity (Anglo): 33.8 oz. Description: 3110-35 Diameter (Anglo): 3.8 in. Inside diameter (anglo) neck: 1.46 in. Inner Diameter (metric) Neck: 37 mm Number Per Carton: 35 Number Per Package: 35 Quantity: 35/pack, 35/carton Type: square Certification/Compliance: USP Class VI, non-cytotoxic, USP It; 661 gt,, pyrogen-free, FDA 21CFR 177.1315 Closure Included: Yes Closure Size: 48 mm Diameter (metric): 98 mm For use with (devices): Long-term storage (application) Graduated: Yes Height (anglo): 7.9 in. Height (metric): 201 mm Opening: Narrow Shape: Square Unit Size: Case of 35

Additional Information

No.	CD-0745
Manufacturer (Brand)	Nalgene
VGKL number	108731103

