



Nalgene® PETG Biotainer Bottles

Product Images



Description

Nalgene PETG Biotainer Bottles Save on cleaning, sterilization and related validation costs compared to reusable bottles. Capacity: 16.9 oz., 500 mL Capacity (metric): 500mL Capacity (anglo): 16.9 oz. Description: 3005-70 Diameter (Anglo): 3 in. Inside Diameter (Anglo) Neck: 1.1 in. Inner Diameter (metric) Neck: 28mm Quantity Per Case: 70 Quantity: 2L HDPE Packer Certified Type: Square Bottle Certification/Compliance: USP Class VI, Non-cytotoxic, USP It;661gt,, Non-pyrogenic, FDA 21CFR 177.1315 Closure Included: Yes Closure Size: 38 mm Diameter (metric): 77 mm For use with (devices): Long-term storage (application) Graduated: Yes Height (anglo): 6.93 in. Height (metric): 176 mm Opening: Narrow Shape: Square Unit Size: Case of 70

Additional Information

No.	CD-0732
Manufacturer (Brand)	Nalgene

