



Nalgene® PVDF True Union Clamps

Product Images



Description

Nalgene PVDF True Union Clamps Use Thermo Scientific Nalgene PVDF union clamps with silicone seals and end caps to repeatedly and reliably seal balloon bottles. Quantity: 1/pk, 1/cs. Type: Tri Clover For use with (devices): Nalgenetrade; sanitary balloon bottles, Nalgenetrade; autoclavable balloon bottles with stopcock Size (Anglo-American) Handle: 3 in. Unit Size: Each

Additional Information

No.	CD-0691
Manufacturer (Brand)	Nalgene

