



## Nalgene® PVDF True Union Clamps

### Product Images

---



16

## Description

Nalgene PVDF True Union Clamps Use Thermo Scientific Nalgene PVDF union clamps with silicone seals and end caps for repeated and reliable sealing of balloon bottles. Quantity: 1/pk, 1/cs. Type: Tri Clover For use with (devices): Nalgenetrade; sanitary balloon bottles, Nalgenetrade; autoclavable balloon bottles with stopcock Size (Anglo-American) Handle: 1.5 in. Unit Size: Each

## Additional Information

|                      |         |
|----------------------|---------|
| No.                  | CD-0690 |
| Manufacturer (Brand) | Nalgene |

